



Placement Guide for the miniSQUAIR®




- Apply parallel to and within 1" of intended incision
- Can be applied to the skin or surgical drapes
- Will adhere to any body contour




Anterior Midline




Posterior Midline




Left Breast



Right Breast




Hip (anterior or lateral incision)




Knee


Placement Guide for the miniSQUAIR®




Right Subcostal




Thoraco-Abdominal




Mid-Sternotomy




Thoracotomy



**Suprapubic
(Pfannenstiel) Incision**



**Perineal
(2 Positions)**



**Perianal
(2 Positions)**