

# miniSQUAIR®

## Accessory Products

### New Product – The "Schultzie"

T-piece connector allows for two smoke capture devices to be attached to a single smoke evacuator.



#### Description:

- For open surgical cases only
- All ports are 22 mm, O.D.; the "arms" connect to 22 mm. I.D. spiral adapters
- The "stem" with the spiral adapter already in place connects to the ULPA filter
- Approved through the "510 (k) non-submission rationale" process
- Patent and trademark applied for



The "Schultzie" allows 2 smoke capture devices to be attached to a single smoke evacuator.

#### Why Perioperative Nurses Want The "Schultzie"

- Easy, quick application of the spiral adapter to the ULPA filter
- No need for two smoke evacuators when two surgeons operate simultaneously with two smoke capture devices (2 "pencils," 2 "mini's" or 1 of each)
- Disposable along with the smoke capture devices

Does Not  
Contain Latex



## Non-Sterile Extension Tubing with Prefilter

72" of tubing connects to the miniSQUAIR® to extend the range for placement of the smoke evacuator.



### Description:

- 25 mm I.D. x 72" corrugated cuffed tubing
- 25 mm x 22 mm reduction fitting
- Pre-filter
- 22 mm. O.D. straight adapter
- 22 mm I.D. x 2" spiral adapter

Approved through the "510 (k) non-submission rationale" process.

Does Not  
Contain Latex



### Ordering Information

Description: T-piece for dual smoke capture devices or The "Schultzie"  
Part Number: SQNS 10011-01  
Packaging: 10 Units/Box, Non-Sterile

Description: **miniSQUAIR Non-Sterile Extension Tubing**  
25mm Inner Diameter x 6 Ft with pre-filter and adapters  
Part Number: SQNS 20018-01  
Packaging: 15 Units/Box, Non-Sterile

## mini SQUAIR®

For a trial of the miniSQUAIR®,  
Extension Tubing with Prefilter, or  
The "Schultzie" call 952-345-1112  
e-mail [info@nascentsurgical.com](mailto:info@nascentsurgical.com)  
or visit [www.nascentsurgical.com](http://www.nascentsurgical.com)

NASCENT  SURGICAL, LLC®

Nascent Surgical, LLC  
6595 Edenvale Blvd., Suite #140  
Eden Prairie, MN 55346